

# LOW TEMPERATURE PLASMA STERILIZER

Sterilization Solution of Plasmapp

**STERLINK**  
PREMIUM LINE



\*plasmapp

“Trusted safeguard for all medical  
devices around the world”

 plasmapp

[www.plasmapp.com](http://www.plasmapp.com)

# Core Technologies for obtaining SAFE and Bio-compatible surface

Plasma is defined as a quasi-neutral gas having a collective behavior (or simply as the 4th state of matter). The plasma discharge voltage is determined by the Paschen's law in which pressure is a key parameter to determine the discharge.



## Under Independently Controllable Pressure

Gentle plasma can be discharged under a vacuum condition without using additional gas.



## Unique Vacuum Technology

Plasma is discharged inside the package by using dielectric barrier of the vacuum state.



## By Discharging Gentle and Powerful Plasma

The pressure can be spatially controlled to discharge gentle and powerful plasma.



## Operated by Validated Reliable Process

Plasma discharge and overall process are validated for obtaining outstanding performance.

## Global Networks

It exports technologies in the world and provides customized solutions that meet the local clients' diversified needs through the global network.

**FDA**  
K212200/K212193/K212198

**CE**  
0068

**TGA**  
323926

**HC**  
102958

**ANVISA**  
81046479001

**SFDA**  
MDMA 19090377

**MFDS**  
18-4069/19-4339

**TFDA**  
DHA05603302400

**TFDA**  
KOR 6209038

**RZN**  
No P3H 2022/17932

**MDA**  
GC5817320

**MHLW**  
JET-302AKBZX

**MDSAP**  
MEDICAL DEVICE SINGLE AUDIT PROGRAM

**TÜV Rheinland**

**GMP**  
GOOD MANUFACTURING PRACTICE

**UL LISTED**  
LABORATORY EQUIPMENT E92914

**RoHS**

# Low temperature plasma sterilizer for wide range of medical devices

The world's first direct injection pouch-type  
**patented sterilization technology** has been applied.  
**SAL of  $10^{-6}$  is validated** to ensure sterilization stability of medical devices.

(SAL: Sterility Assurance Level, 99.9999%)

**STERLINK** *FPS-15s plus*



**STERLINK** *mini*





Low Temperature



Easy Operation



Fast Sterilization



Compact Design




Lumen Sterilization



Real-Time Monitoring

**STERLINK** *lab*







## Plasma sterilization solution is applicable to wide range of medical specialties.

Dental

Plastic Surgery

Ophthalmology

Dermatology

Veterinary

Orthopedics

ENT

### Various Diagnosis & Surgeries

With the development of technology, the surgical instruments and diagnostic equipment have been more diversified.

### Inefficient Sterilization Methods



The medical instruments got complicated and minimized following the trend. However, autoclave can't sterilize heat/moisture sensitive medical devices. EtO Gas Sterilizer requires high maintenance and ethylene oxide has potential hazard.

### Efficient Instruments Manage

It is important to increase the efficiency to deal with patients quickly. Increasing medical instruments turnover rate could be one of the best ways to increase the efficiency

# Novel sterile packaging to enable super-fast 7 min. sterilization cycle

\* STERPACK



- Direct sterilant injection enhances sterilization efficiency.
- Direct vacuum pumping shortens purification process time.
- Independent pressure control allows unique sterilization cycle including compression process.
- Compression maximizes sterilization efficiency to obtain 7 minutes cycle.
- Impermeable pouch allows direct sterilization with compression process and vacuum sealing after the cycle.
- Vacuum sealing visualizes sterile condition which extends shelf life of sterilized devices and provides advanced infection control solution.



Super-fast sterilization cycle



Low temperature reducing maintenance cost of expensive device (Operating temperature < 55°C)



15 FDAs certificates including CE



Reliable sterilization performance (SAL:  $10^{-6}$  with compliance ISO 14937)

## Sterilization Performance

STERLINK verifies sterilization performance through the following lumen tests

### Single-channel lumen claims for STERLINK

- $\varnothing$  0.7 x 500 mm Stainless steel
- $\varnothing$  2.0 x 1,500 mm Stainless steel
- $\varnothing$  1.0 x 2,000 mm PTFE

STERLINK – inside diameter 1mm and a length up to **2,000mm**

Competitor A  
length up to 500mm

Competitor B  
length up to 1,000mm

Plasmapp  
**STERLINK**







# Innovative technologies with patents and research paper



## Forced convective heating for low-temperature sterilization

Publication expected in 2021 (AIP) - <https://doi.org/10.1063/5.0048688>


Review of Scientific Instruments 92, 064902 (2021) conducted by Dr. Youbong Lim, Dr. Wonho Choe, Dr. Seung Hun Lee, Dr. Jun Young Kim and Dr. Hyun Jeong Jeon

### Conclusions

A novel sterile pouch using an impermeable film is presented to improve the heating process, and the experimental and numerical investigations are performed to find that the convective heat transfer coefficient for the forced convective heating is increased more than five times when compared to the natural convective heating. By virtue of the impermeability, the sterilization process is also improved to obtain an overall sterilization cycle completed within 7.5 min, and it is interpreted in terms of the pouch compression.

## Instrument Tracking System (ITS)

**Equipment tracking using barcode**  
Real-time monitoring to check the STERLINK



## Comparison Chart

Sterilizer	Temperature	Cycle Time	Sterilant	Feature
Autoclave	Up to 134°C	60 min + 1 hour cooldown	Hot steam	<ul style="list-style-type: none"> <li>Cloth sterilization</li> <li>Risk of burn</li> <li>Long cycle time</li> <li>Damage of medical instruments</li> </ul>
E.O. Gas	Up to 60°C	Over 120 min + 12 hour ventilation	E.O. gas	<ul style="list-style-type: none"> <li>Highly dangerous toxic gas</li> <li>Long cycle time</li> </ul>
Plasma Sterilizer	Up to 60°C	70 min	H <sub>2</sub> O <sub>2</sub>	<ul style="list-style-type: none"> <li>High cost</li> <li>Additional ventilation</li> <li>Large volume</li> <li>Very low efficiency</li> </ul>
STERLINK	Up to 60°C	Pouch Mode : 7min Pouch Plus Mode : 14 min Chamber Mode : 36 min	H <sub>2</sub> O <sub>2</sub> direct injection	<ul style="list-style-type: none"> <li>Fast sterilization cycle</li> <li>Easy maintenance</li> <li>Economic cost</li> <li>ITS™ system</li> <li>Compact size</li> <li>Eco-friendly</li> <li>Ergonomic design</li> <li>Safeness of heat sensitive instrument</li> </ul>



# Your most effective surgeries with STERLINK sterilizer system

The first non-US plasma sterilizer clearance by US FDA



510(k) clearance for STERLINK sterilizer system of Plasmapp (K212200 / K212193 / K212198)



## Low temperature plasma chamber sterilization Various sizes of medical devices can be sterilized.

FPS-15s plus

mini

lab

14ℓ


264 x 410 x 125

7ℓ

206 x 346 x 113


58ℓ

330 x 700 x 250



# Compatible Medical Devices

- Cranial pressure transducer cables
- Cryoprobes
- Defibrillator paddles
- Dopplers
- Electrocautery instruments
- Endoscopic instruments
- Esophageal dilators
- Fiberoptic lite cables
- Laryngoscope blades
- Laser handpieces, fibers, and accessories
- Metal instruments
- Ophthalmic lenses (diagnostic, magnifying)
- Patient lead cables
- Pigmentation handpieces
- Radiation therapy equipment
- Resectoscope / working elements and sheaths
- Rigid endoscopes
- Shaver handpieces
- Flexible endoscopes
- Stereotactic equipment and batteries
- Trocar sheaths
- Ultrasound probes
- Video cameras and couplers
- Micro instrument kit
- Laparoscopic Grasping Forceps
- Bipolar Forceps with cables
- Vessel Sealer
- Electrosurgical Pencils
- Endotracheal Tube
- Surgical Staples
- Micro Instrument Set
- Oxygen tubing
- Suction & Irrigation tube
- Oxygen mask
- Surgical tips
- Catheters
- Biopsy forceps





# Compatible Materials

Metal	<ul style="list-style-type: none"> <li>• Aluminum</li> <li>• Stainless Steel</li> </ul>
Non-metal	<ul style="list-style-type: none"> <li>• Acrylonitrile butadiene styrene (ABS)</li> <li>• Delrin</li> <li>• Ethyl vinyl acetate (EVA)</li> <li>• Fluorinated ethylene propylene (PFTE)</li> <li>• Glass, USP Type 1 borosilicate</li> <li>• High density polyethylene</li> <li>• Kraton</li> <li>• Latex</li> <li>• Polyethylene</li> <li>• Monel</li> <li>• Phenolic resin</li> <li>• Polycarbonate (PC)</li> <li>• Polyetherimide (PEI)</li> <li>• Polyethylene terephthalate (PET)</li> <li>• Polymethyl methacrylate (PMMA)</li> <li>• Polyphenyl sulphone</li> <li>• Polypropylene (PP)</li> <li>• Polystyrene</li> <li>• Polysulfone</li> <li>• Polytetrafluoroethylene (PTFE)</li> <li>• Polyurethane</li> <li>• Polyvinyl chloride (PVC)</li> <li>• Silicon</li> </ul>

## Non-Compatible Materials

### Absorptive Materials (Cellulose)

- Wood
- Linen
- Paper
- Sponge



### Liquids (Moisture)



### Oil

### Powders






# STERLINK *FPS-15s plus*



Pouch Mode



Chamber Mode



STERPACK & STERPACK plus



STERLOAD



## Specification

Size (W x D x H)	433 x 614 x 437 mm	Vacuum Pump	Pump built-in type
Chamber (W x D x H)	264 x 410 x 125 mm (14 ℓ)	Weight	67kg
Diagonal Length (Chamber)	470 mm	Sterilant Cassette	1 Cycle

Mode	Cycle time	Capacity	Sterilant
Pouch Mode	STERPACK: 7 min	1 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.1 ml/cell)
Pouch Plus Mode	STERPACK Plus: 14 min	4.5 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.3 ml/cell)
Chamber Mode	STERLOAD: 36 min	14 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.9 ml/cell)



# STERLINK *mini*



Pouch Mode



Chamber Mode



STERPACK




STERLOAD mini



## Specification

Size (W x D x H)	275 x 440 x 330 mm	Vacuum Pump	Pump stand-alone type
Chamber (W x D x H)	206 x 346 x 113 mm (7 ℓ)	Weight	20kg (Pump module : 21kg)
Diagonal Length (Chamber)	360 mm	Sterilant Cassette	1 Cycle

Mode	Cycle time	Capacity	Sterilant
Pouch Mode	STERPACK: 7 min	1 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.1 ml/cell)
Advanced Mode	STERLOAD mini: 18 min	7 ℓ	H <sub>2</sub> O <sub>2</sub> 58%—59.5% (0.7 ml/cell)



# STERLINK<sup>lab</sup>



Chamber Mode



UV & Dry Mode



STERLOAD lab




58ℓ

## Specification

Size (W x D x H)	780 x 845 x 485 mm	Vacuum Pump	Pump built-in type
Chamber (W x D x H)	330 x 700 x 250 mm (58 ℓ)	Weight	126 kg
Diagonal Length (Chamber)	770 mm	Sterilant Cassette	1 Cycle

Mode	Cycle time	Capacity	Sterilant
Dry Mode	12 min	58 ℓ	-
Sterilization Mode	39 min		H <sub>2</sub> O <sub>2</sub> 58%—59.5% (5 ml/cell)





# Sterilant & Consumables and Accessories

for Medical Devices Sterilization Performance

## Sterilant cassettes

※ 1 cycle use only



**STERPACK**  
• Size : 135 x 280 mm  
• 50 ea / Box

※ 1 cycle use only




**STERPACK plus**  
• Size : 250 x 430 mm  
• 30 ea / Box

※ 1 cycle use only



**STERLOAD lab**  
• 20 ea / Box

※ 1 cycle use only



**STERLOAD**  
• 30ea / Box

※ 1 cycle use only



**STERLOAD mini**  
• 30ea / Box

## Consumables

※ 1 cycle use only



**Tyvek® 100**  
• Size : 100 x 400 mm  
• 120 ea / Box

※ 1 cycle use only



**Tyvek® 200**  
• Size : 200 x 400 mm  
• 90 ea / Box

※ 1 cycle use only



**Tyvek® 300**  
• Size : 300 x 400 mm  
• 60 ea / Box


※ 1 cycle use only



**Sterilization Wrap**  
• Size : 60 x 60 cm  
• 500 ea / Box

※ Tyvek® is a Dupont™ registered trademark.


## Tray




**STERPACK Tray**  
• Size : 195 x 80 x 30 mm



**STERPACK Plus Tray**  
• Size : 260 x 160 x 50 mm



**Chamber Tray (7ℓ)**  
• Size : 289 x 183 x 90.5 mm



**Chamber Tray (14ℓ)**  
• Size : 327 x 252 x 112 mm

## Accessories



### STERSEAL

- Size : 505 x 250 x 145 mm
- Weight : 12 kg
- Temperature : 135°C / 150°C



### Cart

- plus / mini / lab



### Label Printer

- Size : 120x 102 x 146 mm
- Weight : 0.5 kg



### Label Sticker Roll

- Label size : 40 x 57 mm

※ 1 cycle use only



### CI tape

- Width : 20 mm
- Length : 50 m

※ 1 cycle use only



### CI strip

- Size : 18 x 105 mm
- 250 ea / Box

※ 1 cycle use only



### BI

- Time : 30 min
- 50 ea / Box

## Maintenance

※ 1 cycle use only



### VC Cartridge

- 1ea
- For the vacuum calibration

※ 1 cycle use only



### HEPA Filter

- 1ea
- Every 6 months

※ 1 cycle use only



### Ozone Filter

- 1ea
- Every 6 months

※ 1 cycle use only



### Chamber Filter

- 1ea
- Every 6 months

※ 1 cycle use only



### FPS-15s Plus Deodorizing Filter

- 1ea
- Every 12 months

※ 1 cycle use only



### mini Deodorizing Filter

- 1ea
- Every 12 months



### Pump Oil (LB100)

- 1ea / 280ml
- Every 12 months



### Oil Replacement Kit

- for FPS-15s Plus & mini
- Every 12 months




### Filter Replacement Kit

- for FPS-15s Plus & mini
- Every 6 months



### FPS-15s Plus Filter Kit

- for FPS-15s Plus
- Every 12 months



### mini Filter Kit

- for mini
- Every 12 months

\*plasmapp

[www.plasmapp.com](http://www.plasmapp.com)

# STERLINK

## PREMIUM LINE

Low Temperature Plasma Sterilizer



**Plasmapp Co., Ltd.**

**Sales** | 13F, 476, Gangnam-daero, Gangnam-gu, Seoul, Republic of Korea [06120]

**R&D 1** | 9, Giheungdanji-ro 24beon-gil, Giheung-gu, Yongin-si, Gyeonggi-do, Republic of Korea [17086]

**R&D 2** | 372, Dongbu-daero, Osan-si, Gyeonggi-do, Republic of Korea [18151]

**Factory** | 102, Cheombok-ro, Dong-gu, Daegu, Republic of Korea [41061]

[www.plasmapp.com](http://www.plasmapp.com) | [sales@plasmapp.com](mailto:sales@plasmapp.com)

Specifications are subject to change without any notice.

Some products listed in this catalog are not available for sales in all markets due to pending regulatory clearance.

Please contact the local sales representative for current products availability.

Copyright© 2017 Plasmapp Co., Ltd. All rights reserved.

**STERPR-2309-EN (Rev.0)**